


TLR4 Mouse anti-Human Monoclonal (PE) (HTA125) Antibody - LS-C179342 - LSBio

CatalogID:	LS-C179342
Target:	toll-like receptor 4 (TLR4)
Synonyms:	TLR4 Antibody, CD284 Antibody, CD284 antigen Antibody, Homolog of Drosophila toll Antibody, TLR-4 Antibody, TOLL Antibody, ARMD10 Antibody, Toll-like receptor 4 Antibody, HToll Antibody
Family / Subfamily:	Toll-like Receptor / not assigned-Toll-like Receptor
Host	TLR4 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	HTA125
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	TLR4 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	TLR4 antibody was raised against human TLR4 transfected cells.
Specificity:	Human TLR4
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 1% BSA, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	20 µg
Concentration:	20 µg/ml

Flow Cytometry Image:

Flow Cytometry



Ba/F3 cells co-transfected
with human TLR4 and MD-2

— Anti-TLR4-PE(D077-5)
- - Isotypic control

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences