

| TLR4/MD-2 Complex Rat anti-Mouse Monoclonal (FITC) (MTS510) Antibody - LS-C179344 - LSBio   |  |
|---|--|
| <b>CatalogID:</b>   | LS-C179344   |
| <b>Target:</b>  | Toll-like Receptor 4 (TLR4) / Lymphocyte antigen 96 (MD-2) Complex   |
| <b>Host</b>   | Rat  |
| <b>Clonality:</b>   | Monoclonal   |
| <b>Isotype:</b>   | IgG2a  |
| <b>Clone Name:</b>  | MTS510   |
| <b>Conjugations:</b>  | Fluorescein (FITC)   |
| <b>Immunogen Species:</b>   | Mouse  |
| <b>Antigen Type:</b>  | Cells  |
| <b>Immunogen:</b>   | Mouse TLR4-MD-2 complex expressing Ba/F3 cells.                      |
| <b>Specificity:</b>   | TLR4/MD-2 Complex  |
| <b>Reactivity:</b>  | Mouse  |
| <b>Purification:</b>  | Protein G purified   |
| <b>Presentation:</b>  | PBS, 1% BSA, 0.09% sodium azide                                      |
| <b>Recommended Storage:</b>   | Store at 4°C.  |
| <b>Usage Summary:</b>   | Flo: 20 ul/test;   |
| <b>Uses:</b>  | Flow Cytometry (Optimal dilution to be determined by the researcher) |
| <b>Size:</b>  | 1,000 µl   |
| <b>Requested From:</b>  | Japan  |
| <b>Laboratory Reagent For In Vitro Research Use Only</b><br>Not for resale without prior written consent from LifeSpan BioSciences, Inc.<br>Created on 9/26/2014<br>© 2014 LifeSpan BioSciences |  |