


TLR4/MD-2 Complex Rat anti-Mouse Monoclonal (PE) (MTS510) Antibody - LS-C179345 - LSBio	
CatalogID:	LS-C179345
Target:	Toll-like Receptor 4 (TLR4) / Lymphocyte antigen 96 (MD-2) Complex
Host	Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	MTS510
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Mouse
Antigen Type:	Cells
Immunogen:	Mouse TLR4-MD-2 transfected Ba/F3 cells.
Specificity:	TLR4/MD-2 Complex
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 1% BSA, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	20 µg
Concentration:	20 µg/ml

Flow Cytometry Image:

Flow Cytometry



Ba/F3 cells co-transfected with mouse TLR4 and MD-2

■ Anti-Mouse TLR4-MD-2-PE(D079-5)

□ Isotypic control

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences