

RP105 / CD180 Mouse anti-Human Monoclonal (PE) (MHR73) Antibody - LS-C179347 - LSBio

CatalogID:	LS-C179347
Target:	CD180 molecule
Synonyms:	CD180 Antibody, CD180 antigen Antibody, Lymphocyte antigen 64 Antibody, Ly78 Antibody, LY64 Antibody, CD180 molecule Antibody, RP105 Antibody
Host	CD180 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	MHR73
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	RP105 / CD180 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	RP105 / CD180 antibody was raised against m12HRP mouse B cell line expressing human RP105.
Specificity:	Human RP105 / CD180
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 1% BSA, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Usage Summary:	Flo: 20 ul/test;
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences