

CD151 Mouse anti-Human Monoclonal (PE) (SFA1.2B4) Antibody - LS-C179350 - LSBio	
CatalogID:	LS-C179350
Target:	CD151 molecule (Raph blood group)
Synonyms:	CD151 Antibody, CD151 antigen Antibody, gp27 Antibody, Membrane glycoprotein SFA-1 Antibody, PETA3 Antibody, RAPH Antibody, SFA-1 Antibody, TSPAN24 Antibody, PETA-3 Antibody, SFA1 Antibody, MER2 Antibody, Tetraspanin-24 Antibody, Tspan-24 Antibody
Family / Subfamily:	Tetraspan / not assigned-Tetraspan
Host	CD151 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	SFA1.2B4
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD151 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD151 antibody was raised against human SFA-1 transfectant.
Specificity:	Human CD151
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 1% BSA, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Usage Summary:	Flo: 20 ul/test;
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	