

CD44 Mouse anti-Human Monoclonal (FITC) (15C6) Antibody - LS-C179403 - LSBio

CatalogID:	LS-C179403
Target:	CD44 molecule (Indian blood group)
Synonyms:	CD44 Antibody, CDW44 Antibody, CD44 antigen Antibody, Cell surface glycoprotein CD44 Antibody, CSPG8 Antibody, ECMR-III Antibody, Epican Antibody, HUTCH-I Antibody, LHR Antibody, MDU2 Antibody, MC56 Antibody, MIC4 Antibody, Hermes antigen Antibody, Hyaluronate receptor Antibody, Phagocytic glycoprotein 1 Antibody, PGP-1 Antibody, PGP-I Antibody, Pgp1 Antibody, Phagocytic glycoprotein I Antibody, MDU3 Antibody, CD44R Antibody, HCELL Antibody, Heparan sulfate proteoglycan Antibody
Host	CD44 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	15C6
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD44 antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD44 antibody was raised against mML-1 human leukemia cells.
Specificity:	human CD44
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 1% BSA, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Uses:	Flow Cytometry (10 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	1,000 µl
Concentration:	50 µg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences