

DNAM-1 / CD226 Rat anti-Human Monoclonal (FITC) (TX25) Antibody - LS-C179430 - LSBio	
<b>CatalogID:</b>	LS-C179430
<b>Target:</b>	CD226 molecule
<b>Synonyms:</b>	CD226 Antibody, Adhesion glycoprotein Antibody, CD226 molecule Antibody, CD226 antigen Antibody, DNAM-1 Antibody, DNAM1 Antibody, DNAX accessory molecule-1 Antibody, TLiSA1 Antibody, PTA1 Antibody, DNAX accessory molecule 1 Antibody
<b>Host</b>	CD226 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	TX25
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	DNAM-1 / CD226 antibody was raised against Human
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	DNAM-1 / CD226 antibody was raised against human Recombinant CD226/DNAM-1.
<b>Specificity:</b>	Human DNAM-1 / CD226
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, 1% BSA, 0.09% sodium azide
<b>Recommended Storage:</b>	Store at 4°C.
<b>Uses:</b>	Flow Cytometry (10 µg/ml) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	50 µg
<b>Concentration:</b>	0.5 mg/ml
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
<b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b>	
Created on 9/26/2014 © 2014 LifeSpan BioSciences	