


DLK1 / Pref-1 Rat anti-Mouse Monoclonal (FITC) (24-11) Antibody - LS-C179443 - LSBio	
CatalogID:	LS-C179443
Target:	delta-like 1 homolog (Drosophila) (DLK1)
Synonyms:	DLK1 Antibody, Delta-like protein dlk Antibody, DLK-1 Antibody, Delta1 Antibody, FA1 Antibody, Dlk Antibody, Preadipocyte factor 1 Antibody, Pref-1 Antibody, Protein delta homolog 1 Antibody, PG2 Antibody, ZOG Antibody, Fetal antigen 1 Antibody, PREF1 Antibody, Secredeltin Antibody
Host	DLK1 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	24-11
Conjugations:	Fluorescein (FITC)
Immunogen Species:	DLK1 / Pref-1 antibody was raised against Mouse
Antigen Type:	Fusion protein
Immunogen:	DLK1 / Pref-1 antibody was raised against pref-1-Fc.
Specificity:	Mouse DLK1 / Pref-1
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 1% BSA, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Uses:	Flow Cytometry (50 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	0.5 mg/ml

Flow Cytometry Image:

■ Flow cytometry



Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences