

BAFF Receptor / BAFF-R Mouse anti-Human Monoclonal (PE) (8A7) Antibody - LS-C179457 - LSBio

CatalogID:	LS-C179457
Target:	tumor necrosis factor receptor superfamily, member 13C (TNFRSF13C)
Synonyms:	TNFRSF13C Antibody, BAFFR Antibody, BR3 Antibody, BAFF receptor Antibody, BLyS receptor 3 Antibody, CD268 Antibody, CD268 antigen Antibody, CVID4 Antibody, Prolixin Antibody, BAFF-R Antibody, BROMIX Antibody
Family / Subfamily:	TNF Receptor / Not assigned-TNF Receptor
Host	TNFRSF13C antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	8A7
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	BAFF Receptor / BAFF-R antibody was raised against Human
Antigen Type:	Cells
Immunogen:	BAFF Receptor / BAFF-R antibody was raised against human CD268/BAFF-R/BR3 transfectant.
Specificity:	Human BAFF Receptor / BAFF-R
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 1% BSA, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences