

**CD52 Rat anti-Mouse Monoclonal (PE) (BTG-2G) Antibody - LS-C179461 - LSBio**

<b>CatalogID:</b>	LS-C179461
<b>Target:</b>	CD52 molecule
<b>Synonyms:</b>	CD52 Antibody, Cambridge pathology 1 antigen Antibody, CD52 antigen Antibody, CD52 molecule Antibody, CDW52 Antibody, CAMPATH-1 antigen Antibody, He5 Antibody
<b>Host</b>	CD52 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	BTG-2G
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Immunogen Species:</b>	CD52 antibody was raised against Mouse
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	CD52 antibody was raised against mouse IL-2R alpha-transgenic mice splenocyte.
<b>Specificity:</b>	Mouse CD52
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, 1% BSA, 0.09% sodium azide
<b>Recommended Storage:</b>	Store at 4°C.
<b>Usage Summary:</b>	Flo: 20 ul/test;
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	50 tst
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences