

TLR4/MD-2 Complex Mouse anti-Mouse Monoclonal (PE) (UT15) Antibody - LS-C179465 - LSBio	
CatalogID:	LS-C179465
Target:	Toll-like Receptor 4 (TLR4) / Lymphocyte antigen 96 (MD-2) Complex
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	UT15
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Mouse
Antigen Type:	Cells
Immunogen:	Mouse TLR4/MD-2 transfected Ba/F3 cells.
Specificity:	TLR4/MD-2 Complex
Reactivity:	Mouse
Purification:	Protein A purified
Presentation:	PBS, 1% BSA, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Usage Summary:	Flo: 20 ul/test;
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	