

FOXP3 Rat anti-Mouse Monoclonal (FITC) (MF-14) Antibody - LS-C179486 - LSBio

CatalogID:	LS-C179486
Target:	forkhead box P3 (FOXP3)
Synonyms:	FOXP3 Antibody, AIID Antibody, DIETER Antibody, Forkhead box protein P3 Antibody, JM2 Antibody, Forkhead box P3 Antibody, FOXP3delta7 Antibody, XPID Antibody, PIDX Antibody, SCURFIN Antibody
Host	FOXP3 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	MF-14
Conjugations:	Fluorescein (FITC)
Immunogen Species:	FOXP3 antibody was raised against Mouse
Antigen Type:	Recombinant protein
Immunogen:	FOXP3 antibody was raised against recombinant Mouse full length Foxp3.
Specificity:	Mouse FOXP3
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 1% BSA, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Usage Summary:	Flo: 20 ul/test;
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences