


**CD45RO Mouse anti-Human Monoclonal (UCHL1) Antibody - LS-C179538 - LSBio**

<b>CatalogID:</b>	LS-C179538
<b>Target:</b>	CD45RO
<b>Host</b>	CD45RO antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a,k
<b>Clone Name:</b>	UCHL1
<b>Immunogen Species:</b>	CD45RO antibody was raised against Human
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	CD45RO antibody was raised against iL-2 dependent T cell line (CA1).
<b>Specificity:</b>	CD45RO
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 50% glycerol
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	IHC: Heat treatment is necessary for paraffin-embedded sections;IP: 2 ug/200-300 ul of cell extract from 5x10E6 cells;
<b>Uses:</b>	IHC - Paraffin (1 - 5 µg/ml), Western blot (10 µg/ml), Immunoprecipitation, Flow Cytometry (1 - 10 µg/ml) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Concentration:</b>	1 mg/ml

**Flow Cytometry Image:**

**Flow Cytometry**



<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	