

CD45RO Mouse anti-Human Monoclonal (FITC) (UCHL1) Antibody - LS-C179539 - LSBio

CatalogID:	LS-C179539
Target:	CD45RO
Host	CD45RO antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	UCHL1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CD45RO antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD45RO antibody was raised against iL-2 dependent T cell line (CA1).
Specificity:	CD45RO
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 1% BSA, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Uses:	Flow Cytometry (10 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	1,000 µl
Concentration:	50 µg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences