

PDGFRB / PDGFR Beta Mouse anti-Human Monoclonal (PE) (18A2) Antibody - LS-C179678 - LSBio

CatalogID:	LS-C179678
Target:	platelet-derived growth factor receptor, beta polypeptide (PDGFRB)
Synonyms:	PDGFRB Antibody, CD140B Antibody, Beta pdgfr Antibody, Beta type PDGFR Antibody, BPDGFR Antibody, CD140b antigen Antibody, IBGC4 Antibody, JTK12 Antibody, PDGF-beta R Antibody, PDGF-R-beta Antibody, PDGFR Antibody, PDGFRbeta Antibody, PDGF receptor beta Antibody, PDGF-R Antibody, PDGFbetaR Antibody, PDGFR Beta Antibody, PDGFR-1 Antibody, PDGFR-beta Antibody, BetaPDGFR Antibody, Pdgf beta-receptor Antibody, PDGFR1 Antibody, PDGFRB/CEV14 fusion gene Antibody
Family / Subfamily:	Protein Kinase / PDGF Receptor
Host	PDGFRB antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	18A2
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	PDGFRB / PDGFR Beta antibody was raised against Human
Antigen Type:	Transfected Cells
Immunogen:	PDGFRB / PDGFR Beta antibody was raised against human CD140b transfected NIH/3T3.
Specificity:	Human PDGFRB / PDGFR Beta
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 1% BSA, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Usage Summary:	Flo: 20 ul/test;
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences