

HLA-A2 Mouse anti-Human Monoclonal (PE) (BB7.2) Antibody - LS-C179715 - LSBio

CatalogID:	LS-C179715
Target:	MHC Class I HLA-A2
Host	MHC Class I HLA-A2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b,k
Clone Name:	BB7.2
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	HLA-A2 antibody was raised against Human
Antigen Type:	Purified protein
Immunogen:	HLA-A2 antibody was raised against papain solubilized HLA-A2.
Specificity:	MHC Class I HLA-A2
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 1% BSA, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Usage Summary:	Flo: 20 ul/test;
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences