


MYC / c-Myc Mouse Monoclonal (PL14) Antibody - LS-C179776 - LSBio	
<b>CatalogID:</b>	LS-C179776
<b>Target:</b>	v-myc avian myelocytomatosis viral oncogene homolog
<b>Synonyms:</b>	MYC Antibody, BHLHe39 Antibody, C-Myc Antibody, MRTL Antibody, Myc proto-oncogene protein Antibody, Proto-oncogene c-Myc Antibody, Rats1 Antibody, Transcription factor p64 Antibody, v-myc Antibody
<b>Host</b>	MYC antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	PL14
<b>Antigen Type:</b>	Fusion protein
<b>Immunogen:</b>	MYC / c-Myc antibody was raised against 6myc-tag fusion protein.
<b>Specificity:</b>	c-Myc
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 50% glycerol
<b>Recommended Storage:</b>	Store at -20°C.
<b>Uses:</b>	ICC (2 µg/ml), Western blot (1 µg/ml), Immunoprecipitation (5 µg/tst) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µg
<b>Concentration:</b>	2 mg/ml

Western Blot Image:

■ Western blotting



Lane 1: Myc-Tagged-p38  
Lane 2: Myc-Tagged-JNK1  
Lane 3: Myc-Tagged-MKK3

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences