


MYC / c-Myc Mouse Monoclonal (Biotin) (PL14) Antibody - LS-C179777 - LSBio

CatalogID:	LS-C179777
Target:	v-myc avian myelocytomatosis viral oncogene homolog
Synonyms:	MYC Antibody, BHLHe39 Antibody, C-Myc Antibody, MRTL Antibody, Myc proto-oncogene protein Antibody, Proto-oncogene c-Myc Antibody, Rats1 Antibody, Transcription factor p64 Antibody, v-myc Antibody
Host	MYC antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	PL14
Conjugations:	Biotin
Antigen Type:	Fusion protein
Immunogen:	MYC / c-Myc antibody was raised against 6myc-tag fusion protein.
Specificity:	c-Myc
Purification:	Protein A purified
Presentation:	Sodium phosphate buffer, pH 8.0, 0.3 M sodium chloride, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Uses:	ICC (1:100), Western blot (1:1000) (Optimal dilution to be determined by the researcher)
Size:	100 µl

Western Blot Image:

■ Western blotting



Lane 1: Myc-Tag heparanase
Lane 2: Myc-Tag cdc25

Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences