


MYC / c-Myc Mouse Monoclonal (HRP) (PL14) Antibody - LS-C179778 - LSBio

CatalogID:	LS-C179778
Target:	v-myc avian myelocytomatosis viral oncogene homolog
Synonyms:	MYC Antibody, BHLHe39 Antibody, C-Myc Antibody, MRTL Antibody, Myc proto-oncogene protein Antibody, Proto-oncogene c-Myc Antibody, Rats1 Antibody, Transcription factor p64 Antibody, v-myc Antibody
Host	MYC antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	PL14
Conjugations:	Horseradish Peroxidase (HRP)
Antigen Type:	Fusion protein
Immunogen:	MYC / c-Myc antibody was raised against 6myc-tag fusion protein.
Specificity:	c-Myc
Purification:	Protein A purified
Presentation:	PBS, Preservative, Stabilizer
Recommended Storage:	Store at 4°C.
Uses:	Western blot (1:1000 - 1:4000) (Optimal dilution to be determined by the researcher)
Size:	100 µl

Western Blot Image:

免疫沈降後のウェスタンブロットングでの使用例
間接法とHRP-DirectTの比較



Requested From:

Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences