

Phosphotyrosine Mouse Monoclonal (FITC) (6D12) Antibody - LS-C179874 - LSBio

CatalogID:	LS-C179874
Target:	Phosphotyrosine
Host	Phosphotyrosine antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	6D12
Conjugations:	Fluorescein (FITC)
Antigen Type:	Hapten
Immunogen:	Phosphotyrosine antibody was raised against kLH conjugated phosphotyrosine.
Specificity:	Phosphotyrosine
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.5 M sodium chloride, 50% glycerol
Recommended Storage:	Store at 4°C.
Uses:	ICC (10 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	100 µl
Concentration:	1 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences