


**BrdU Mouse Monoclonal (2B1) Antibody - LS-C179876 - LSBio**

<b>CatalogID:</b>	LS-C179876
<b>Target:</b>	Bromodeoxyuridine (BrdU)
<b>Host</b>	Bromodeoxyuridine (BrdU) antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	2B1
<b>Antigen Type:</b>	Hapten
<b>Immunogen:</b>	BrdU antibody was raised against 5-Iodouridine-Ovalbumin.
<b>Specificity:</b>	BrdU
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 50% glycerol
<b>Recommended Storage:</b>	Store at -20°C.
<b>Uses:</b>	IHC (10 µg/ml), ICC, Flow Cytometry (10 µg/ml) (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Concentration:</b>	1 mg/ml

**Immunohistochemistry Image:**

**■ Immunohistochemistry**



Mouse testis

Paraffin section

**Requested From:** Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences