


BrdU Mouse Monoclonal (PE) (2B1) Antibody - LS-C179877 - LSBio	
CatalogID:	LS-C179877
Target:	Bromodeoxyuridine (BrdU)
Host	Bromodeoxyuridine (BrdU) antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	2B1
Conjugations:	Phycoerythrin (PE)
Antigen Type:	Hapten
Immunogen:	BrdU antibody was raised against 5-Iodouridine-Ovalbumin.
Specificity:	BrdU
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, 1% BSA, 0.09% sodium azide
Recommended Storage:	Store at 4°C.
Usage Summary:	Flo: 20 ul/test;
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	50 tst
<p>Flow Cytometry Image:</p> <p>■ Flow cytometry</p> <p>Bromodeoxyuridine incorporated in nuclei in Raji cells</p>  <p> Anti-Bromodeoxyuridine Isoypic control </p>	
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences