

**GALP / Alarin Rabbit anti-Human Polyclonal (Cy3) Antibody - LS-C182045 - LSBio**

<b>CatalogID:</b>	LS-C182045
<b>Target:</b>	galanin-like peptide (GALP)
<b>Synonyms:</b>	GALP Antibody, Alarin Antibody, Gal-like peptide Antibody, Galanin-like peptide Antibody
<b>Host</b>	GALP antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Cy3
<b>Immunogen Species:</b>	GALP / Alarin antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	GALP / Alarin antibody was raised against synthetic peptide corresponding to human Alarin.
<b>Specificity:</b>	Human Alarin.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified
<b>Presentation:</b>	Liquid
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences