

IL22RA1 / IL22R Rat anti-Mouse Monoclonal (aa18-228) (PerCP) (12K118) Antibody - LS-C182892
- LSBio

CatalogID:	LS-C182892
Target:	interleukin 22 receptor, alpha 1 (IL22RA1)
Synonyms:	IL22RA1 Antibody, CRF2-9 Antibody, IL-22RA1 Antibody, IL-22 receptor subunit alpha-1 Antibody, IL22R Antibody, IL-22R-alpha-1 Antibody, ZcytoR11 Antibody, IL22R1 Antibody, Interleukin 22 receptor Antibody
Family / Subfamily:	Interleukin
Host	IL22RA1 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	12K118
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	IL22RA1 / IL22R antibody was raised against Mouse
Antigen Type:	Recombinant protein
Immunogen:	IL22RA1 / IL22R antibody was raised against recombinant corresponding to aa18-228 from mouse IL22RA1 expressed in CHO cell line.
Specificity:	Mouse IL22RA1.
Epitope:	aa18-228
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	Saline, 0.5% BSA, 0.1% sodium azide
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences