

PVRL3 / Nectin-3 Rat anti-Mouse Monoclonal (aa50-400) (FITC) (13A62) Antibody - LS-C183001 - LSBio	
CatalogID:	LS-C183001
Target:	poliovirus receptor-related 3 (PVRL3)
Synonyms:	PVRL3 Antibody, CDW113 Antibody, DKFZP566B0846 Antibody, Nectin 3 Antibody, PPR3 Antibody, Poliovirus receptor-related 3 Antibody, CD113 Antibody, CD113 antigen Antibody, NECTIN-3 Antibody, PVRR3 Antibody
Host	PVRL3 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	13A62
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PVRL3 / Nectin-3 antibody was raised against Mouse
Immunogen:	PVRL3 / Nectin-3 antibody was raised against extracellular domain (aa50-400) from mouse Nectin-3 fused to IgG Fc
Specificity:	Mouse Nectin-3. Does not cross-react with human or rat.
Epitope:	aa50-400
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, 0.1% BSA, 0.02% sodium azide
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Usage Summary:	Immunohistochemistry: Frozen (PFA or acetone/methanol fixation)
Uses:	IHC - Frozen, ICC, Immunoprecipitation, Flow Cytometry, ELISA, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	100 µg
Concentration:	0.1 mg/ml
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	