

FGF21 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C183182 - LSBio	
CatalogID:	LS-C183182
Target:	fibroblast growth factor 21 (FGF21)
Synonyms:	FGF21 Antibody, Fibroblast growth factor 21 Antibody, FGF-21 Antibody
Family / Subfamily:	HBGF / not assigned-HBGF
Host	FGF21 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	FGF21 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	FGF21 antibody was raised against recombinant corresponding to human Fibroblast growth factor 21.
Specificity:	Human FGF21.
Reactivity:	Human
Purification:	Ammonium sulfate precipitation
Presentation:	10 mM PBS, PH 7.4, 0.03% Proclin 300, 50% glycerol
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	50 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	