

NeuN Rabbit anti-Mouse Polyclonal (Cy3) Antibody - LS-C183399 - LSBio

CatalogID:	LS-C183399
Target:	Neuron Specific Nuclear Protein (NeuN)
Host	Neuron Specific Nuclear Protein (NeuN) antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Cy3
Immunogen Species:	NeuN antibody was raised against Mouse
Antigen Type:	Recombinant protein
Immunogen:	NeuN antibody was raised against recombinant corresponding to mouse NeuN fused to GST-tag.
Specificity:	Mouse NeuN.
Reactivity:	Mouse, Rat
Purification:	Affinity purified
Presentation:	PBS, 0.1% sodium azide, 15 mg/ml BSA
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ICC (1:100) (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences