

ANXA1 / Annexin A1 Rabbit anti-Human Polyclonal (aa2-26) (Cy3) Antibody - LS-C183752 - LSBio	
CatalogID:	LS-C183752
Target:	annexin A1 (ANXA1)
Synonyms:	ANXA1 Antibody, Annexin A1 Antibody, Annexin I Antibody, Annexin-1 Antibody, ANX1 Antibody, Calpactin II Antibody, Chromobindin-9 Antibody, Lipocortin I Antibody, Annexin I (lipocortin I) Antibody, Calpactin-2 Antibody, LPC1 Antibody
Host	ANXA1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Cy3
Immunogen Species:	ANXA1 / Annexin A1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ANXA1 / Annexin A1 antibody was raised against synthetic peptide corresponding to aa2-26 of human Ac2-26.
Specificity:	Human Ac2-26.
Epitope:	aa2-26
Reactivity:	Human
Purification:	Purified
Reconstitution:	Sterile distilled water. Possible additional volumes for resuspension: 100 µl
Presentation:	Lyophilized
Recommended Storage:	Lyophilized powder may be stored at -20°C. Stable for 12 months at -20°C. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Reconstituted product is stable for 12 months at -20°C.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	