

Endothelin Rabbit anti-Human Polyclonal (aa109-125) (FITC) Antibody - LS-C183846 - LSBio	
<b>CatalogID:</b>	LS-C183846
<b>Target:</b>	Endothelin
<b>Host</b>	Endothelin antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	Endothelin antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	Endothelin antibody was raised against synthetic peptide corresponding to aa109-125 of human Endothelin-1.
<b>Specificity:</b>	Human Endothelin-1.
<b>Epitope:</b>	aa109-125
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified
<b>Presentation:</b>	Liquid
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µl
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	