

GHRL / Ghrelin-Obestatin Rabbit anti-Human Polyclonal (aa52-75) (Rhodamine) Antibody - LS-C183870 - LSBio

CatalogID:	LS-C183870
Target:	Ghrelin/Obestatin Prepropeptide (GHRL)
Synonyms:	GHRL Antibody, Ghrelin Antibody, Growth hormone secretagogue Antibody, Motilin-related peptide Antibody, Obestatin Antibody, Appetite-regulating hormone Antibody, MTLRP Antibody, Protein M46 Antibody
Host	GHRL antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Rhodamine
Immunogen Species:	GHRL / Ghrelin-Obestatin antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	GHRL / Ghrelin-Obestatin antibody was raised against synthetic peptide corresponding to aa52-75 of human Ghrelin.
Specificity:	Human Ghrelin.
Epitope:	aa52-75
Reactivity:	Human
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences