

GHRL / Ghrelin-Obestatin Rabbit anti-Rat Polyclonal (Cy3) Antibody - LS-C183875 - LSBio	
CatalogID:	LS-C183875
Target:	Ghrelin/Obestatin Prepropeptide (GHRL)
Synonyms:	GHRL Antibody, Ghrelin Antibody, Growth hormone secretagogue Antibody, Motilin-related peptide Antibody, Obestatin Antibody, Appetite-regulating hormone Antibody, MTLRP Antibody, Protein M46 Antibody
Host	GHRL antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Cy3
Immunogen Species:	GHRL / Ghrelin-Obestatin antibody was raised against Rat
Antigen Type:	Synthetic peptide
Immunogen:	GHRL / Ghrelin-Obestatin antibody was raised against synthetic peptide corresponding to rat Ghrelin.
Specificity:	Rat Ghrelin.
Reactivity:	Rat, Mouse
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	