

**GLP2 Rabbit anti-Human Polyclonal (aa1-34) (Cy3) Antibody - LS-C183888 - LSBio**

<b>CatalogID:</b>	LS-C183888
<b>Target:</b>	Glucagon-like Peptide 2 (GLP2)
<b>Host</b>	Glucagon-like Peptide 2 (GLP2) antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Cy3
<b>Immunogen Species:</b>	GLP2 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	GLP2 antibody was raised against synthetic peptide corresponding to aa1-34 of human GLP-2.
<b>Specificity:</b>	Human GLP-2.
<b>Epitope:</b>	aa1-34
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified
<b>Presentation:</b>	Liquid
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences