

| MCHR2 Rabbit anti-Human Polyclonal (aa237-251) (Rhodamine) Antibody - LS-C183962 - LSBio | |
|------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| CatalogID: | LS-C183962 |
| Target: | melanin-concentrating hormone receptor 2 (MCHR2) |
| Synonyms: | MCHR2 Antibody, G protein-coupled receptor slt Antibody, GPRv17 Antibody, GPR145 Antibody, G protein-coupled receptor 145 Antibody, G-protein coupled receptor 145 Antibody, MCH receptor 2 Antibody, MCHR-2 Antibody, MCH-2R Antibody, MCH-R2 Antibody, MCH2 Antibody, Mch2 receptor Antibody, SLT receptor Antibody, MCH2R Antibody, SLT Antibody |
| Family / Subfamily: | GPCR / Melanin-concentrating hormone receptor |
| Host | MCHR2 antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Rhodamine |
| Immunogen Species: | MCHR2 antibody was raised against Human |
| Antigen Type: | Synthetic peptide |
| Immunogen: | MCHR2 antibody was raised against synthetic peptide corresponding to aa237-251 of human MCH-2R. |
| Specificity: | Human MCH-2R. |
| Epitope: | aa237-251 |
| Reactivity: | Human |
| Purification: | Purified |
| Presentation: | Liquid |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light. |
| Uses: | ELISA (Optimal dilution to be determined by the researcher) |
| Size: | 100 µl |
| Requested From: | Japan |
| Laboratory Reagent For In Vitro Research Use Only | |
| Not for resale without prior written consent from LifeSpan BioSciences, Inc. | |
| Created on 9/26/2014 © 2014 LifeSpan BioSciences | |