

MC4R / Melanocortin 4 Receptor Rabbit anti-Human Polyclonal (aa306-332) (FITC) Antibody - LS-C183965 - LSBio	
CatalogID:	LS-C183965
Target:	melanocortin 4 receptor (MC4R)
Synonyms:	MC4R Antibody, Melanocortin receptor 4 Antibody, MC4-R Antibody, MC4 Antibody, MC4 receptor Antibody, Melanocortin 4 receptor Antibody
Family / Subfamily:	GPCR / Melanocortin
Host	MC4R antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	MC4R / Melanocortin 4 Receptor antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MC4R / Melanocortin 4 Receptor antibody was raised against synthetic peptide corresponding to aa306-332 of human MC4 Receptor.
Specificity:	Human MC4 Receptor.
Epitope:	aa306-332
Reactivity:	Human
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	