

MC4R / Melanocortin 4 Receptor Rabbit anti-Rat Polyclonal (aa306-332) (FITC) Antibody - LS-C183968 - LSBio

CatalogID:	LS-C183968
Target:	melanocortin 4 receptor (MC4R)
Synonyms:	MC4R Antibody, Melanocortin receptor 4 Antibody, MC4-R Antibody, MC4 Antibody, MC4 receptor Antibody, Melanocortin 4 receptor Antibody
Family / Subfamily:	GPCR / Melanocortin
Host	MC4R antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	MC4R / Melanocortin 4 Receptor antibody was raised against Rat
Antigen Type:	Synthetic peptide
Immunogen:	MC4R / Melanocortin 4 Receptor antibody was raised against synthetic peptide corresponding to aa306-332 of rat MC4 Receptor.
Specificity:	Rat MC4 Receptor.
Epitope:	aa306-332
Reactivity:	Rat
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences