

ITLN1 / Omentin Rabbit anti-Human Polyclonal (aa163-199) (FITC) Antibody - LS-C184019 - LSBio	
<b>CatalogID:</b>	LS-C184019
<b>Target:</b>	intelectin 1 (galactofuranose binding) (ITLN1)
<b>Synonyms:</b>	ITLN1 Antibody, Endothelial lectin HL-1 Antibody, Galactofuranose-binding lectin Antibody, HL1 Antibody, INTL Antibody, Intelectin-1 Antibody, ITLN-1 Antibody, LFR Antibody, HL-1 Antibody, Omentin Antibody, HIntL Antibody, ITLN Antibody
<b>Host</b>	ITLN1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	ITLN1 / Omentin antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	ITLN1 / Omentin antibody was raised against synthetic peptide corresponding to aa163-199 of human Intelectin.
<b>Specificity:</b>	Human Intelectin.
<b>Epitope:</b>	aa163-199
<b>Reactivity:</b>	Human
<b>Purification:</b>	Purified
<b>Presentation:</b>	Liquid
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	