

Orexin A Rabbit anti-Human Polyclonal (aa16-33) (Cy3) Antibody - LS-C184025 - LSBio

CatalogID:	LS-C184025
Target:	Orexin A
Host	Orexin A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Cy3
Immunogen Species:	Orexin A antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	Orexin A antibody was raised against synthetic peptide corresponding to aa16-33 of human Orexin A.
Specificity:	Human Orexin A.
Epitope:	aa16-33
Reactivity:	Human, Mouse, Rat, Pig, Sheep
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences