

OSTN / Musclin / Osteocrin Rabbit anti-Rat Polyclonal (aa31-77) (Rhodamine) Antibody - LS-C184037 - LSBio

CatalogID:	LS-C184037
Target:	osteocrin (OSTN)
Synonyms:	OSTN Antibody, MUSCLIN Antibody, Osteocrin Antibody
Host	OSTN antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Rhodamine
Immunogen Species:	OSTN / Musclin / Osteocrin antibody was raised against Rat
Antigen Type:	Synthetic peptide
Immunogen:	OSTN / Musclin / Osteocrin antibody was raised against synthetic peptide corresponding to aa31-77 of rat Musclin.
Specificity:	Rat Musclin.
Epitope:	aa31-77
Reactivity:	Rat
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences