

OSTN / Musclin / Osteocrin Rabbit anti-Human Polyclonal (aa83-132) (Rhodamine) Antibody - LS-C184046 - LSBio	
CatalogID:	LS-C184046
Target:	osteocrin (OSTN)
Synonyms:	OSTN Antibody, MUSCLIN Antibody, Osteocrin Antibody
Host	OSTN antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Rhodamine
Immunogen Species:	OSTN / Musclin / Osteocrin antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	OSTN / Musclin / Osteocrin antibody was raised against synthetic peptide corresponding to aa83-132 of human Musclin.
Specificity:	Human Musclin.
Epitope:	aa83-132
Reactivity:	Human
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	