

PACAP Rabbit anti-Human Polyclonal (aa6-27) (Cy5) Antibody - LS-C184051 - LSBio

CatalogID:	LS-C184051
Target:	adenylate cyclase activating polypeptide 1 (pituitary) (ADCYAP1)
Synonyms:	ADCYAP1 Antibody, PACAP Antibody
Family / Subfamily:	Ligand / not assigned-Ligand
Host	ADCYAP1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Cy5
Immunogen Species:	PACAP antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PACAP antibody was raised against synthetic peptide corresponding to aa6-27 of human PACAP.
Specificity:	Human PACAP.
Epitope:	aa6-27
Reactivity:	Human, Rat, Sheep
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences