

PROK1 / EG-VEGF Rabbit anti-Human Polyclonal (aa81-105) (Cy3) Antibody - LS-C184092 - LSBio	
CatalogID:	LS-C184092
Target:	prokineticin 1 (PROK1)
Synonyms:	PROK1 Antibody, EG-VEGF Antibody, EGVEGF Antibody, Mambakine Antibody, PRK1 Antibody, Prokineticin-1 Antibody, PK1 Antibody, Prokineticin 1 Antibody
Host	PROK1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Cy3
Immunogen Species:	PROK1 / EG-VEGF antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PROK1 / EG-VEGF antibody was raised against synthetic peptide corresponding to aa81-105 of human Prokineticin 1.
Specificity:	Human Prokineticin 1.
Epitope:	aa81-105
Reactivity:	Human
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	