

SORL1 Rabbit anti-Human Polyclonal (aa2177-2214) (FITC) Antibody - LS-C184128 - LSBio	
CatalogID:	LS-C184128
Target:	sortilin-related receptor, L(DLR class) A repeats containing (SORL1)
Synonyms:	SORL1 Antibody, C11orf32 Antibody, LR11 Antibody, LRP9 Antibody, gp250 Antibody, Mosaic protein LR11 Antibody, SorLA-1 Antibody, SORLA Antibody, Sortilin-related receptor Antibody
Host	SORL1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	SORL1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	SORL1 antibody was raised against synthetic peptide corresponding to aa2177-2214 of human SORL1.
Specificity:	Human SORL1.
Epitope:	aa2177-2214
Reactivity:	Human
Purification:	Purified
Presentation:	Liquid
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	