

TRPM8 Rabbit anti-Rat Polyclonal (aa656-680) Antibody - LS-C184151 - LSBio

CatalogID:	LS-C184151
Target:	transient receptor potential cation channel, subfamily M, member 8 (TRPM8)
Synonyms:	TRPM8 Antibody, CMR1 Antibody, Cold-menthol receptor type 1 Antibody, LTrpC-6 Antibody, TRPP8 Antibody, Trp-p8 Antibody, LTRPC6 Antibody
Family / Subfamily:	Ion Channel / Calcium channel - TRP
Host	TRPM8 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Immunogen Species:	TRPM8 antibody was raised against Rat
Antigen Type:	Synthetic peptide
Immunogen:	TRPM8 antibody was raised against synthetic peptide corresponding to aa656-680 of rat TRPM8.
Specificity:	Reconizes rat TRPM8.
Epitope:	aa656-680
Reactivity:	Rat
Purification:	Purified
Reconstitution:	Sterile distilled water. Possible additional volumes for resuspension: 200 µl
Presentation:	Lyophilized
Recommended Storage:	Lyophilized powder may be stored at -20°C. Stable for 12 months at -20°C. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Reconstituted product is stable for 12 months at -20°C.
Uses:	ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences