

PECAM-1 / CD31 Mouse anti-Human Monoclonal (PerCP) (11n31) Antibody - LS-C184192 - LSBio	
CatalogID:	LS-C184192
Target:	platelet/endothelial cell adhesion molecule 1 (PECAM1)
Synonyms:	PECAM1 Antibody, CD31 Antibody, CD31 antigen Antibody, EndoCAM Antibody, PECA1 Antibody, PECAM-1 Antibody, GPIIA Antibody
Host	PECAM1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	11n31
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	PECAM-1 / CD31 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	PECAM-1 / CD31 antibody was raised against recombinant corresponding to CD31/PECAM-1.
Specificity:	Human CD31.
Reactivity:	Human
Purification:	Protein G purified
Presentation:	1ml saline, 0.5% BSA, 0.1% sodium azide
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	