

| CXCR6 Rat anti-Mouse Monoclonal (aa1-351) (PerCP) (221002) Antibody - LS-C184213 - LSBio | |
|--|---|
| CatalogID: | LS-C184213 |
| Target: | chemokine (C-X-C motif) receptor 6 (CXCR6) |
| Synonyms: | CXCR6 Antibody, CD186 antigen Antibody, CXCR-6 Antibody, STRL33 Antibody, TYMSTR Antibody, BONZO Antibody, CD186 Antibody, CDw186 Antibody, CXCR6 Antibody |
| Family / Subfamily: | GPCR / Chemokine |
| Host | CXCR6 antibody was produced in Rat |
| Clonality: | Monoclonal |
| Isotype: | IgG2b |
| Clone Name: | 221002 |
| Conjugations: | Peridinin-chlorophyll-protein (PerCP) |
| Immunogen Species: | CXCR6 antibody was raised against Mouse |
| Antigen Type: | Recombinant protein |
| Immunogen: | CXCR6 antibody was raised against recombinant corresponding to aa1-351 from mouse CXCR6 expressed in CHO cells (Q9EQ16). |
| Specificity: | Mouse CXCR6. |
| Epitope: | aa1-351 |
| Reactivity: | Mouse |
| Purification: | Protein A/G purified |
| Presentation: | Supplied in saline, 0.5% BSA, 0.09% sodium azide |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 12 months. Sensitive to light. |
| Uses: | Flow Cytometry (Optimal dilution to be determined by the researcher) |
| Size: | 100 tst |
| Requested From: | Japan |
| Laboratory Reagent For In Vitro Research Use Only | |
| Not for resale without prior written consent from LifeSpan BioSciences, Inc. | |
| Created on 9/26/2014 | |
| © 2014 LifeSpan BioSciences | |