

OVOL1 Mouse anti-Human Monoclonal (460C1a) Antibody - LS-C184697 - LSBio

CatalogID:	LS-C184697
Target:	ovo-like zinc finger 1 (OVOL1)
Synonyms:	OVOL1 Antibody, HOVO1 Antibody, Ovo-like 1(Drosophila) Antibody, Ovo homolog-like 1 Antibody
Host	OVOL1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	460C1a
Immunogen Species:	OVOL1 antibody was raised against Human
Immunogen:	OVOL1 antibody was raised against recombinant Human Ovo-Like 1(Drosophila).
Specificity:	Human OVOL1
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 1% BSA, 0.05% sodium azide
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Uses:	Dot Blot (Optimal dilution to be determined by the researcher)
Size:	50 µg
Concentration:	100 µg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences