

CO-029 / TSPAN8 Rat anti-Mouse Monoclonal (PE) (11n186) Antibody - LS-C184905 - LSBio	
CatalogID:	LS-C184905
Target:	tetraspanin 8 (TSPAN8)
Synonyms:	TSPAN8 Antibody, CO-029 Antibody, TM4SF3 Antibody, Tspan-8 Antibody, Tetraspanin 8 Antibody, Tetraspanin-8 Antibody
Family / Subfamily:	Tetraspan / not assigned-Tetraspan
Host	TSPAN8 antibody was produced in Rat
Clonality:	Monoclonal
Clone Name:	11n186
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CO-029 / TSPAN8 antibody was raised against Mouse
Antigen Type:	Recombinant protein
Immunogen:	CO-029 / TSPAN8 antibody was raised against recombinant corresponding to mouse TSPAN8 expressed in CHO cells.
Specificity:	Mouse TSPAN8.
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	1ml saline, 0.5% BSA, 0.1% sodium azide
Recommended Storage:	May be stored at 4°C. For long-term storage, aliquot and store at 4°C. Do not freeze. Aliquots are stable for at least 12 months. Sensitive to light.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	