

TNFRSF17 / BCMA Rat anti-Human Monoclonal (FITC) (Vicky-1) Antibody - LS-C18662 - LSBio	
<b>CatalogID:</b>	LS-C18662
<b>Target:</b>	tumor necrosis factor receptor superfamily, member 17 (TNFRSF17)
<b>Synonyms:</b>	TNFRSF17 Antibody, B cell maturation antigen Antibody, BCMA Antibody, CD269 antigen Antibody, B-cell maturation factor Antibody, B-cell maturation protein Antibody, BCM Antibody, CD269 Antibody
<b>Family / Subfamily:</b>	TNF Receptor / not assigned-TNF Receptor
<b>Host</b>	TNFRSF17 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	Vicky-1
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	TNFRSF17 / BCMA antibody was raised against Human
<b>Immunogen:</b>	TNFRSF17 / BCMA antibody was raised against recombinant human BCMA extracellular domain (aa 2-54)
<b>Specificity:</b>	Recognizes human BCMA. Does not cross-react with mouse BCMA.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, 0.02 % sodium azide.
<b>Recommended Storage:</b>	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light. Store undiluted.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	