

CD80 / B7 Mouse anti-Bovine Monoclonal (RPE) (IL-A159) Antibody - LS-C187501 - LSBio	
CatalogID:	LS-C187501
Target:	CD80 molecule
Synonyms:	CD80 Antibody, Activation B7-1 antigen Antibody, B7-1 Antibody, B7.1 Antibody, BB1 Antibody, CD28LG Antibody, CD28LG1 Antibody, CD80 antigen Antibody, LAB7 Antibody, CD80 molecule Antibody, CTLA-4 counter-receptor B7.1 Antibody
Host	CD80 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	IL-A159
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	CD80 / B7 antibody was raised against Bovine
Reactivity:	Bovine, Sheep
Purification:	Protein G purified
Reconstitution:	Distilled water
Presentation:	PBS, 0.09% sodium azide, 1% BSA, 5% sucrose
Recommended Storage:	Store at +4°C. Do not freeze. Protect from light.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/26/2014	
© 2014 LifeSpan BioSciences	