

CD86 Mouse anti-Bovine Monoclonal (RPE) (IL-A190) Antibody - LS-C187504 - LSBio

CatalogID:	LS-C187504
Target:	CD86 molecule
Synonyms:	CD86 Antibody, Activation B7-2 antigen Antibody, B7-2 Antibody, B7.2 Antibody, B70 Antibody, CD86 molecule Antibody, CTLA-4 counter-receptor B7.2 Antibody, CD86 antigen Antibody, FUN-1 Antibody, LAB72 Antibody, BU63 Antibody, CD28LG2 Antibody, Ctla-4 counter-receptors b7.2 Antibody
Host	CD86 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	IL-A190
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	CD86 antibody was raised against Bovine
Reactivity:	Bovine, Sheep
Purification:	Protein G purified
Reconstitution:	1 ml sterile distilled water
Presentation:	PBS, 0.09% sodium azide, 1% BSA, 5% sucrose
Recommended Storage:	Store at +4°C. Do not freeze. Protect from light.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/26/2014

© 2014 LifeSpan BioSciences